






	Dew 350	S2	S1	9 Plus	Hand Plus
					
MODEL #	MM010599	MM011305	MM011091	MM011343	MM010964
WARRANTY	2-year motor 90-day accessories	2-year motor 90-day accessories	2-year motor 90-day accessories	1-year motor 90-day accessories	90-day
AVERAGE VACUUM STRENGTH	0-275mmHg (+/- 10%)	50-280mmHg (+/-15%)	50-280mmHg (+/- 15%)	0-260mmHg (+/- 15%)	0-300mmHg
2-PHASE MODE MASSAGE EXPRESSION	No	Yes	Yes	Yes	No
CYCLE LEVEL	One-Dial Button to Adjust suction strength	5-level (38, 42, 46, 50, 54)	5-level (38, 42, 46, 50, 54)	One touch adjusts suction strength and cycle speed	N/A
VACUUM LEVEL	N/A	Massage L1-L5 Expression L1-L12	Massage L1-L5 Expression L1-L12	Massage L1-L5 Expression L1-L10	N/A
FLANGE SIZE AVAILABLE	20mm, 24mm, 28mm, 32mm	20mm, 24mm, 28mm, 32mm	20mm, 24mm, 28mm, 32mm	20mm, 24mm, 28mm, 32mm	24mm only
BASIC ACCESSORY KIT	Pump Unit W/ Attached Power Cord 24mm Breast Flanges (2) Wide Neck Bottles (2) Bottle Locking Rings(2) Bottle Discs (2) Bottle Cap (2) Duckbill Valve (2) Filter (4) Filter Tubing(2) Tubing (2)	S2 Plus Motor Unit (1) Power Cord/AC Adapter (1) 24mm Breast Flanges (2) 28mm Breast Flange (2) Wide Neck Bottles (2) Bottle Locking Rings (2) Bottle Discs (2) Bottle Caps (2) Bottle Locking Rings (2) Bottle Discs (2) Bottle Caps (2) Duckbill Valves (2) Backflow Protectors (2) Tubings (2)	S1 Plus Motor Unit (1) Power Cord/AC Adapter (1) 24mm Breast Flanges (2) 28mm Breast Flange (2) Wide Neck Bottles (2) Bottle Locking Rings (2) Bottle Discs (2) Bottle Caps (2) Duckbill Valves (2) Backflow Protectors (2) Tubings (2)	9 Plus Motor Unit (1) Power Cord/AC Adapter (1) 24mm Breast Flanges (2) Wide Neck Bottles (2) Bottle Locking Rings (2) Bottle Discs (2) Bottle Caps (2) Duckbill Valves (2) Backflow Protectors (2) Tubings (2)	24mm flange (1) Narrow neck bottle(1) Tubing(1) Membrane(1) Duckbill valve (1)
CLOSED SYSTEM	Yes	Yes	Yes	Yes	Yes
LCD DISPLAY	No	Yes	Yes	Yes	No
NIGHT LIGHT	No	2 levels	2 levels	No	No
POWER	AC Electrical outlet	AC Electrical outlet	AC/Electrical outlet Built-in Rechargeable Battery	AC/Electrical outlet Built-in Rechargeable Battery	Manual
INPUT OUTPUT POWER	220V / 50-60Hz	100V - 240V / 50-60 Hz 12V, 2A	100V - 240V / 50-60 Hz 12V, 2A	100V - 240V / 50-60 Hz 12V, 2A	N/A
WEIGHT	4 lbs	2.5 lbs	3 lbs	0.5 lb	0.3 lb
COLOR	Pink	Pink	Blue	White	Clear
DIMENSION	8" X 5.5"	7.5" x 7"	7.5" x 7"	4.4"X 3" X 2"	N/A
WIRELESS USAGE (BUILT-IN BATTERY)	N/A	N/A	3 hrs	2 hrs	N/A